



PerkinElmer Introduces TurboMass™ 5.4 GC/MS Software

PerkinElmer introduces new TurboMass™ 5.4 GC/MS software, now expanded to deliver enhanced productivity, robustness and easier navigation.

TurboMass™ GC/MS software is designed for a wide variety of application needs. It provides features and functionalities such as QA/QC method evaluation, environmental CLP-like reporting, customized reporting and much, much more. It is robust and powerful enough to meet the most demanding laboratory requirements. In addition, it is extremely easy to use and has the flexibility and power of sophisticated data-evaluation tools.

TurboMass™ 5.4 software is available as an upgrade to current users of PerkinElmer TurboMass™ Gold or Clarus® GC/MS systems, or to users looking to purchase a new Clarus® GC/MS. Combining the latest in TurboMass™ GC/MS software with the benefits of the Clarus® 600 GC/MS instruments provides users the finest in bench-top GC/MS technology.

For more information, please visit www.perkinelmer.com/TurboMass